


Pocket Guide for Repositioning Sling

Supine Lifting (Using Repositioning Sling):

1. Apply repositioning sling by log rolling patient.
2. Ensure main (padded) spreader bar is PERPENDICULAR to patients shoulders and located just above the abdomen.
3. On each side of patient, attach two lower sets of loops to spreader bar as shown (use green loops).
4. On each side of patient, attach two upper sets of loops to spreader bar as shown (use yellow loops).
5. Attach set of head loops if needed (red or yellow). Lift patient.



Patient Turning (Using Repositioning Sling):

1. Ensure spreader bar is PERPENDICULAR to patients shoulders and is located just above the abdomen and attach 2 lower sets of loops to the center hook of spreader bar as shown (use red loops).
2. Attach 2 upper sets of loops to the center hook of spreader bar as shown (use red loops).
3. Press up button on hand controller to roll patient.



-WARNINGS-

Always Ensure Bed Is In Locked Position before Using Lift

Visibly inspect sling prior to each use to ensure sling is the correct type, size and design to handle lifting; the sling is not damaged, torn, worn, discolored or past its useful life; that the sling's straps are correctly attached to the spreader bar; and that the sling is tested with resident in it at a few inches over bed prior to actual lifting.

Never leave patient/resident unattended in sling.